
Automated ambulance disinfection

(aHP = aerosolized Hydrogen Peroxide)

Never again dangerous gaps in your surface and final disinfections!!!

Do any of the following situations affect your emergency rescue service or your company in general?

You have a shortage of skilled workers?

You have budget shortages?

You have a lack of time?

You suffer from disinfection frustration?

Are you worried about your vehicle and work
safety?

Then why haven't you implemented
Disinfection Automation in your ambulance
vehicles yet?

What is your EXCUSE?

Dramatic basic problems and worrying hygiene trends

Budget Scarcity: No financial resources for "expensive" sanitation and disinfection works

Time Shortage: No time resources for challenging final disinfections

Personnel shortage: No qualified personnel to carry out ("*shortage of skilled workers*")

Lack of Will: Supervisors and hygiene managers ignore innovative disinfection procedures

**The Hygiene dilemma:
A vicious circle**

Problem view instead of problem-solving joy:
"*Naysayers*" see only problems instead of concrete solutions

Rising infection risks:
Multi-resistant germs, fungi, viruses and spores

Counseling Resistance: New solutions are not heard and maybe not wanted

Human Factor: Lack of focus on automated, validatable disinfection techniques

Nevertheless, so far only the manual scrubbing and wiping method method is used in most final disinfections.

According to "DGKH" (German Society for Hospital Hygiene), this manual, human-based disinfection method achieves an efficiency of approx. 70%. **What about the remaining 30%?**

Conclusion: That is definitely not enough!

Application: Healthcare

Application: Food & pharma production

Meat processing

Lock systems

Food production

Veterinary medicine


Special applications:

Air rescue

- ✓ Sporicide Disinfection of Air Transfer Isolators
- ✓ Maximum protection of employees and patients
- ✓ Preparation via MF Professional System in disinfection chambers / treatment rooms

Step by step application of automated vehicle disinfection

Always keep an eye on the relative humidity values

Simple overdose and adaptation possible!

Example:

The calculated volume of the room to be disinfected amounts to 150 m³.

The value of the rel. humidity is 57%.

Therefore, 25% of the value must be added to the settings for 150m³ of the DiosolGenerator:

$$150\text{m}^3 + 37,5\text{m}^3 \approx \mathbf{190\text{m}^3}$$

190m³ is the setting at which the disinfection is to be conducted for a rel. Humidity of 57%.

Value of the rel. Humidity %	DiosolGenerator setting
<30	+ 100% of the volume value
30-40	+ 50% of the volume value
40-55	Not required
55-80	+ 25% of the volume value
>80	+ 50% of the volume value

Special accessories for comfortable vehicle preparation


Demand-oriented applications:

- ✓ Device in the patient area
- ✓ Device via syringe application between driver and patient area

Investment & running costs in the rescue service

Net price calculation with DiosolGenerator MF Professional & Diosol 19 for rescue vehicles

*Acquisition cost

The **DiosolGenerator** is being applied with the disinfectant Diosol (Hydrogen peroxide with or without silver ions) in vehicles by way of aerosol disinfection („*Nebulization*“ or „*Cold fogging*“).

Producing a fine mist (droplet size between 2-10 μm), the Diosol automatically covers all surfaces, the air and the smallest gaps and spaces.

The disinfection measure

- Generation of a disinfecting micro-aerosol with an even distribution within the vehicle. You may now control the complete vehicle architecture!
- Immediate proof of efficacy by chemical indicators / data loggers on site!
- Wetting of all surfaces down to the farthest fugue
- leaves neither wiped / sticky film nor toxic residues
- low consumption (about 4 ml per m³ vehicle)
- Compatible with sensitive electronics & materials

The Diosol aHP method...

Should not be mixed up with the following technologies:

- Formalin room fumigation (formaldehyde)
- Hydrogen peroxide room fumigation (e.g. Vaporized Hydrogen Peroxide, VHP short: Steris method)
- Hydrogen Peroxide Evaporation (Hydrogen Peroxide Vapor: Bioquell Process)
- The ozone odor neutralization

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Process times & protection


Simple formula: 90 minutes (exposure time) + x minutes (ventilation)

<ul style="list-style-type: none"> • Diosol 12 • Diosol 12 PURE 	12%	 	See Diosol 6	90 minutes	See Diosol 8	<ul style="list-style-type: none"> • With ventilation after 25-30 minutes • Without ventilation 5-7 hours
<ul style="list-style-type: none"> • Diosol 19 • Diosol 19 PURE 	19%	  	See Diosol 6	90 minutes	Additional protection necessary	<ul style="list-style-type: none"> • Without ventilation about 10-12 hours • With adequate ventilation you may re-enter room 35-45 minutes after exposure time

The microbiological effect of the procedure was and is being tested in a multiplicity of reports & studies.

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Diosol disinfectants tested according to EN standards!

EN 1040, EN 1276, EN 1650, EN 14347,
EN 13704, EN 14348, EN 14476, EN 1655,
EN 1657, EN 14204

Depending on country-based regulations for
biocides or medical products!

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

Up to sporicidal efficacy

Microbiological effectiveness of H₂O₂ room disinfection with DiosolGenerators

Advantages of aHP versus „classical“ H₂O₂ gaseous approaches

Effective area ABCD already
by using H₂O₂ < 20%.
Application according to
German BioStoffV

Lower dangerous substance class
(UN 2984 instead of UN 2014)

Much lower PPM load (max.
350 ppm instead of about
1500 ppm)


No surface and material
changes with proper use

High efficacy also against
biofilms

100% electronic compatible (also
being used for disinfection of sensitive
medical equipment)

Saving time, budget, staff and
chemistry. Significantly lower
acquisition costs

**Advantages of
aerosolized
Hydrogen Peroxide
(aHP) vs. Gaseous
systems (VHP, HPV
etc.)**


Special hygiene and instruction
training for new customers.
Permanent customer support by
disinfection specialists

Easy-to-use, sporicidal
validation using "Closed
bioindicators"

Upgrading and further developments
according to customer request

Quicker disinfection method
(Process time, dispersion time
and exposure time)

Electronical compatibility of H₂O₂ nebulization

Time frame for room disinfection

Timeline of an aHP room disinfection with DIOP DiosolGenerators

According to application field, the following options must be chosen:

- DiosolGenerator Model
- Diosol- / Diosol PURE Concentration (3-19%)
- Chemical indicator test sticks (25 vs. 100 ppm)
- In case of validation: DioFog- vs. DioSpore-Controller

Ventilation time

by turning on the room ventilation systems or by venting via open windows or doors. According to Diosol concentration between 0 and 45 minutes

Room preparation time

of room disinfection:
According to field of application about **3-25 minutes**

- measure rel. humidity
- place Generator
- place indicators
- Seal openings with gastight tape
- Cover smoke detector covers
- specific preparation according to field of application
- in case of service disinfections: Putting on protective clothing

Process time

(Nebulization starts, Micro Aerosol formation):
According to device type:

- up to 30 minutes (270 m³)
- up to 60 minutes (450 m³)
- up to 120 minutes (850 m³)

Dispersion time

(Micro Aerosol spreads during and after process time gapless in complete environment):

According to room volume approximately:
up to minute 60 (when using a 270 m³ generator)

Exposure time

Microbiological activity takes **90 minutes** in total including process time and dispersion time (Exception for 850 m³ systems: 120 Minutes exposure time as system works up to 120 minutes)

Follow-up work

Room disinfection:
According to field of application about **3 to 15 minutes**

- Collection of indicators
- Remove gastight seal and smoke detector covers
- specific action according to field of application

Special and optional disinfection equipment for vehicle decontamination

Standard components
(for new national customers
including delivery by car):

Gastight tape

Expert knowledge
certificate for
application

Chemical indicators
(H₂O₂ test sticks)

Diosol canister
of your choice
(3-19%)

Nozzle channel
cleaner

Optional for Multifunctional systems (MF + Tube models):

V4A-
Bottle holder

3 meter
hose including
all connections

Suction hose
with V4A
coupling

1 kg
empty bottle

DIOP
GmbH & Co KG
Disinfection is our Passion

**DIOP Starter
Package for H₂O₂
Room disinfection**
by way of aerosolized
Hydrogen Peroxide
(aHP)
www.disinfection-shop.com

Optional according to
your specific
requirements:

V4A Throttle valve
(e.g. for incubators)

Window foams
for vehicle
disinfection

Universal
smoke detector
cover + telescopic rod

Hygrometer

Bioindicators for
bactericidal
validation (bacteria)


Bioindicators for
sporicidal
validation (spores)

Complicated and insufficient validation options


Agar contact plates


Open germ carriers / spore carriers

Uncomplicated validation options:
DioFog and DioSpore-Controller
bioindicators

Example: Sporicidal validation in emergency vehicles

Thus, an optimal disinfection is guaranteed and a validation can be performed.

A **validation certificate** can be issued by independent reviewers, experts & laboratories.

A validation of the scrubbing and wiping method is not possible.

In many areas it has been recognized that the current hygiene standard is no longer sufficient and institutions as well as companies are going the "new" way of automated vehicle disinfection by using the " H_2O_2 cold fogging" technology.


The DioProtection disinfection process has been used successfully for many years in:

- Pharmaceuticals (laboratories, clean rooms), animal husbandry and food industry
- Medical practices and dental surgeries
- Clinics and hospitals
- Nursing homes
- Rescue and fire departments
- Passenger ships: Crusaders and luxury liners

Hygiene references (small abstract)

Further information in our online shop for hygiene and disinfection

Do you want more exclusive information + a discount offer for flexible vehicle disinfection?

www.hygiene-standard.com

DIOP
GmbH & Co KG
Disinfection is our Passion

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© Springer-Verlag Berlin Heidelberg 2013